

Dr. Odin®
DIGITAL THERMOMETER
Instruction Manual
Model No. : ODT201BT

Intended Use: It is a device used to record body temperature.

Warning :

- ⚠ Read instructions thoroughly before using digital thermometer.
- ⚠ Choking Hazard: Thermometer cap and battery may be fatal if swallowed. Do not allow children to use this device without parental supervision.
- ⚠ Do not use thermometer in ear. Designed use is for oral, rectal, and armpit (axilla) readings only.
- ⚠ Do not place thermometer battery near extreme heat as it may explode.
- ⚠ Note: Use of the probe cover may result in a 0.1°C (0.2°F) discrepancy from actual temperature.
- ⚠ It is recommended the performance should be checked every two years.
- ⚠ Remove battery from the device when not in operation for a long time.
- ⚠ The use of temperature readings for self-diagnosis is dangerous. Consult your doctor for the interpretation of results. Self-diagnosis may lead to the worsening of existing disease conditions.
- ⚠ Do not attempt measurements when the thermometer is wet as inaccurate readings may result.
- ⚠ Do not bite the thermometer. Doing so may lead to breakage and/or injury.
- ⚠ Do not attempt to disassemble or repair the thermometer. Doing so may result in inaccurate readings.
- ⚠ After each use, disinfect the thermometer especially in case the device is used by more than one person.
- ⚠ Do not force the thermometer into the rectum. Stop insertion and abort the measurement when pain is present. Failure to do so may lead to injury.
- ⚠ Do not use thermometer orally after being used rectally.
- ⚠ For children who are two years old or younger, please do not use the devices orally.

PLEASE READ CAREFULLY BEFORE USING


This digital thermometer provides a quick and highly accurate reading of an individual's body temperature. To better understand its functions and to provide years of dependable results, please read all instructions first.

CONTENTS

1 Thermometer, 1 Owners Manual, 1 Storage Case

PRODUCT ILLUSTRATION


Figure 1



PRECAUTION

- * The performance of the device may be degraded should one or more of the following occur:
 - Operation outside the manufacturer's stated temperature and humidity range.
 - Storage outside the manufacturer's stated temperature and humidity range.
 - Mechanical shock (for example, drop test).
 - Patient temperature is below ambient temperature.
- * Portable and mobile RF communications such as microwave ovens and cell phones, which generate strong electrical or electromagnetic fields, near the device, as they may become interference source to the device.

SYMBOL EXPLANATION



SPECIFICATIONS

Type:	Digital Thermometer (Not Predictive)
Range:	32.0°C - 42.9°C (90.0°F -109.9°F) (°C /°F chosen by manufacturer)
Accuracy:	±0.1°C, 35.5°C - 42.0°C (±0.2 °F, 95.9°F -107.6°F) ±0.2°C under 35.5°C or over 42.0°C (±0.4°F under 95.9°F or over 107.6°F) At standard room temperature of 25°C (77.0°F)
Display:	Liquid crystal display, 3 1/2 digits
Memory:	For storing the last measured value
Battery:	One 1.5 V DC, button battery (size LR41)
Battery life:	Approx 200 hours
Dimension:	12.3x1.7x0.9 cm (LxWxH)
Weight:	Approx. 10 grams including battery
Storage Condition	Keep it in cool & dry place.
Alarm (Beeper Mode)	Approx. 10 seconds sound when peak temperature reached

Bluetooth requirements

1. The thermometer requires a bluetooth device with:
 - * Bluetooth 4.0 or later
 - * Android 6.0 or later
 - * IOS 10.0 or later.
2. And works with: . iphone , iPod, iPad . Android Phones and Tablets.


°C / °F SWITCHABLE

Temperature readings are available in the Celsius or Fahrenheit scale (°C/°F; located in the upper right corner of LCD.) With the unit off, press and hold the *On/Off* Button for approximately 3 seconds to change the current setting.

MEMORY MODE

1. Power on the thermometer.
2. During the last memory displays, press and hold the *On/Off* button until entering into Memory Mode.
3. Press the button again to cycle check the last 10 memories.
4. Hold down the on/Off button for 3 seconds or wait for one minute to exit the memory mode.

BATTERY REPLACEMENT

1. Replace battery when “  ” appears in the lower right corner of LCD display.
2. Pull battery cover off as shown in Figure 3.
3. Gently pull out plastic circuit board with battery chamber approximately 1 cm (slightly less than 1/2"). (See Figure 4)
4. Use pointed object such as a pen to remove old battery. Discard battery lawfully. Replace with new 3.0V DC button type CR2032. Be sure battery is installed with “-” polarity facing up. (See Figure 5)
5. Slide battery chamber back into place and attach cover.

Note: The following schematic diagram of battery replacement is from the back of the thermometer




Figure 3



Figure 4



Figure 5

DIRECTIONS

1. Press the *On/Off* Button next to LCD display. A tone will sound as the screen shows *1888.8*, followed by last recorded temperature. After showing the self-test temperature, the thermometer is now in the testing mode.

*Note: If the measured temperature is less than 32 °C or 90 °F, the LCD will display Lo, if the measured temperature is higher than 42.9°C or 109.9°F, the LCD will display Hi.

2. Position thermometer in desired location (mouth, rectum, or armpit.)
 - a) **Oral Use:** Place thermometer under tongue as indicated by “ ✓ ” position shown in Figure 2. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air. Normal temperature between 35.7°C and 37.3°C (96.3 °F and 99.1 °F)
 - b) **Rectal Use:** Lubricate silver probe tip with petroleum jelly for easy insertion. Gently insert sensor approximately 1cm (less than 1/2") into rectum. Normal temperature between 36.2°C and 37.7°C (97.2 °F and 99.9 °F).
 - c) **Armpit Use:** Wipe armpit dry. Place probe in armpit and keep arm pressed firmly at side. From a medical viewpoint, this method will always provide inaccurate readings, and should not be used if precise measurements are required. Normal temperature between 35.2°C and 36.7°C (95.4°F and 98.1°F).
3. The degree sign flashes throughout the testing process. When flashing stops an alarm will beep for approximately 10 seconds. The measured reading will appear on the LCD simultaneously. The minimum measurement time until the signaling tone (beep) must be maintained without exception. The measurement continues even after the buzzer notification. The displayed temperature will not change when thermometer is removed from its testing position.




Figure 2

*Note: Normally the buzzes are “ Bi-Bi-Bi-Bi-” ; Alarm beeps more rapidly when temperature reaches 37.8°C (100 °F) or higher, and the buzzes are “Bi-Bi-Bi----- Bi-Bi-Bi----- Bi-Bi-Bi”

4. To prolong battery life, press the *On/Off* Button to turn unit off after testing is complete. If no action is taken, the unit will automatically shut off after around 10 minutes.
5. Store the thermometer in its protective case.

How to Connect Your Dr. Odin Digital Thermometer via Bluetooth

Step 1: Download the App

First, you'll need to download the official Dr. Odin app. You can find it by searching for "**Dr. Odin**" on the Google Play Store for Android devices or the Apple App Store for iOS devices.

Step 2: Enable Bluetooth

Once the app is installed, open it on your device. The app requires Bluetooth to function, so make sure it's enabled in your smartphone's Settings menu. This allows the app and the Digital Thermometer to communicate with each other.

Step 3: Log In or Create an Account

You'll be prompted to either create a new user account or log in with an existing one. Follow the on-screen instructions to proceed.

Step 4: Pair Your Device

On the app's home screen, navigate to the "**Connect Your Devices**" section and select Digital Thermometer. Next, start your Dr. Odin Digital Thermometer to begin a temperature reading. The app will automatically detect the device. Then tap on device to connect it with your app.

Step 5: Sync Your Data

After the Digital Thermometer completes the measurement, the results will be automatically reflected on the app via Bluetooth. The app will then save this data to your memory for future reference and tracking.

Important Note: To ensure a successful data transfer, keep the app open and active on your phone until the measurement is complete and the results have appeared.

Troubleshooting

If you encounter any issues with the Bluetooth connection, try the following steps:

1. Restart the app.
2. Re-connect the device by repeating the pairing process.

Remember, even if Bluetooth communication fails, you can always view your measurement data directly on the Digital Thermometer's display.

Passim Lifesciences Ltd. warrants this units free of any defect in material, workmanship and operation, under normal use, for a period of one year from the date of purchase. If any part become defective during the warranty period Passim Lifesciences Ltd. will repair or replace (if not repairable) the same, free of cost. This shall not apply to any parts that are considered as expandable or deteriorable in the course of normal use. Passim Lifesciences Ltd. shall be relieved of any liability and warranty shall cease to apply if :

- This is not used in accordance with the instructions in the operational manual.
- It is used with any equipment not complying with the specification of this unit.
- It is not regularly maintained.
- The unit is disassembled, repaired or operated by person not authorized by Passim Lifesciences Ltd.
- Damage caused due to negligency.
- The unit is operated in corrosive materials or in the harmful atmosphere.
- The warranty card is not filled completely and produced at the time of warranty claim.

CLEANING AND DISINFECTION

Wipe the thermometer with a soft clean cloth.

For stubborn stains, wipe the thermometer with a cloth that has been dampened with water or a neutral detergent solution and then wring thoroughly. Finish by wiping with a soft dry cloth.

Observe the following to prevent damage to the thermometer.

- Do not use benzene, thinner, gasoline or other strong solvents to clean the thermometer.
- Do not soak the sensing section in alcohol for long periods of time or attempt to sterilize it using hot water (water at a temperature of 50 °C (122 °F) or higher).
- Do not use ultrasonic washing to clean the thermometer.

Passim Lifesciences Ltd. warrants this units free of any defect in material, workmanship and operation, under normal use, for a period of 6 months from the date of purchase. If any part become defective during the warranty period Passim Lifesciences Ltd. will repair or replace (if not repairable) the same, free of cost. This shall not apply to any parts that are considered as expandable or deteriorable in the course of normal use. Passim Lifesciences Ltd. shall be relieved of any liability and warranty shall cease to apply if :

- This is not used in accordance with the instructions in the operational manual.
- It is used with any equipment not complying with the specification of this unit.
- It is not regularly maintained.
- The unit is disassembled, repaired or operated by person not authorized by Passim Lifesciences Ltd.
- Damage caused due to negligency.
- The unit is operated in corrosive materials or in the harmful atmosphere.
- The warranty card is not filled completely and produced at the time of warranty claim.

Manufactured By:

Passim Lifesciences Ltd., Plot No. 45, Industrial Area, Phase 2, Panchkula, Haryana 134113, India

Mfg. Lic. No.: MFG/MD/2022/000615

WPC ETA Reg. No.: ETA-SD- 20250807420

For complaints or suggestions :

Customer Care No.: 1800 303 3009

(Timing: 9am - 7pm, Mon-Sat.)

Email ID: customercare@droidin.in

Website: www.droidin.in

WARRANTY CARD	
Model No.	Lot No.
Invoice No.	Date of Purchase.....
Purchased By:	Contact No.
Address.....	
Dealer's Name.....	
Dealer's Sign & Stamp	

* Terms & Conditions Apply