



User Manual
CF570BT

Weighing Scale
Body Fat Analyzer



One-click sync with



Apple Health



Google Fit

Mobile&Device



iOS 13.0



Android 6.0

For First Use

1. Battery Installation

Open the battery door, put into batteries with correct polarity.

2. App Installation

Scan the QR code on the above, or download “FitGo” from App Store or Android Google Play.

3. Registration

Turn on bluetooth on mobile device, run the App“FitGo”, follow the App instructions.

Attention

1.Do not use the scale on people who have body implants such as a pacemaker, artificial limbs, contra-
ceptive devices, or metal plates or screws. It may cause the device to malfunction or produce aninaccu-
rate result. When in doubt, consult your physician.

2. Please make sure good contact between bare feet and the metal electrodes, otherwise it will display
your weight only.

3. Ensure the scale is on a flat, even and firm surface. Avoid carpet or soft surfaces.

4. The surface of the scale will be slippery if wet. Please make sure both the glass and your feet are dry
before use.

5. This scale is not water proof. Do not expose to moisture.

6. Clean the scale with a slightly damp cloth. DO not use any chemical/abrasive cleaning agents to clean
the scale.

7. Please remove the batteries, if the scale is not to being used for a long time.

8. The appliance is intended solely for domestic use within the home.

9. Children should be supervised to ensure that they do not play with the appliance.

10. When weighing for the first time, there is a big difference in weight reading.

The electronic scale will be automatically adjusted when it is used for the first time or after being moved.
The weighed data includes the weight of the prod uct itself, and accurate data can be obtained by weigh-
ing again.

11. Why are the weighing data different on different grounds?

Because of the highly sensitive sensor systems of electronic scales, weighing needs to be placed on a flat
and firm surface (such as uneven ground or placed on a carpet or foam floor mat will cause the weighing
data to be inaccurate)

12. Reasons for inconsistent weighing data


Unsteady posture: body shaking, etc. make it difficult to lock the data during weighing. It is recommend-
ed to stand steadily in the center of the electronic scale

Ground factor: If the ground is uneven or placed on an unstable surface such as a carpet, foam pad, etc.,
the weighing data will be incorrect.

An important point- It is recommended to stand steadily in the center of the electronic scale.

The actual features and interface are subject to the latest update.Features are subject to
change without prior notice in future updates.

24 parameters in 1 measurment



Multi Users

8 users share one user account in the device.

Target weight

Users can set their own target weight to customize their own health plan.

Imported & Marketed By:
Passim Medichem Agencies
186, Ind. Area Phase II, Chandigarh - 160 002
Corporate Office:
902, 9th Floor, WallFort House, S.V. Road,
Near Citi Mall, Goregaon (W), Mumbai-400 062
WPC ETA Reg. No.: ETA-SD-20260505026
For any Complaint/ Suggestion
please contact Customer Care No.: 1800 309 3009
(Timing: 9am -7pm, Mon. - Sat.)
Email ID: customercare@drodin.in, **Website:** www.drodin.in